

AutoLumo *A*1860

Compact and intelligent with more accurate results



Up to 180 Tests/Hour

25 Reagents Positions

60 Sample Positions

Single Cuvette (RV), Single Test

Compact but Versatile

- Surrounding washing by vacuum.
- Stepping workstation for magnetic separation.
- Innovative vibration mixing.

Easy operation and save time

- Iconized management of consumables, automatic switch of dual-substrate.
- Keep cooling for reagents after switch off, random access.
- Auto-retesting, RFID reading.

More accurate results

- Lower carryover, lower failure rate.
- Lower background and higher sensitivity.
- 1.8s to mix well with higher precision.

A 1860: L x W x H = 1200 mm x 635 mm x 736 mm, Net Weight = 183 kg

Autobio Diagnostics Co., Ltd.

No.87, Jingbei Yi Rd, National Eco & Tech Zone, Zheng zhou City, China-450016

Facebook: @autobioGlobal, Youtube Channel: Autobio Diagnostics,

E-mail: info@autobio-diagnostics.com, Whatsapp: +8618595767835, Tel: +86-371-6798-5313



Tumor Markers

AFP
CEA
tPSA
CA50
fPSA
CA125
CA15-3
CA19-9
Ferritin
β 2-Microglobulin
NSE
CYFRA 21-1
SCCA
CA72-4
Pepsinogen I
Pepsinogen II
CA242
HE4
ProGRP
PIVKA-II
G-17*
HER-2/neu*
S100*

Vasculitis

Anti-PR3 IgG*
Anti-GBM IgG*
Anti-MPO IgG*

Infetious Diseases

Anti-HCV
HAV-IgM
HEV-IgM
HEV-IgG
Anti-HIV
HIV Ag/Ab Combo
Anti-TP
HP-IgG
Anti-HTLV*
VCA-IgA*
EA-IgA*
VCA-IgM*
NA1-IgA*
EA-IgG*

Respiratory Disease

Influenza A IgM
Influenza B IgM
MP-IgM
MP-IgG
CP-IgM
LP-IgM
RSV-IgM
ADV-IgM
COXB-IgM
EV71-IgM

Bone Metabolism

Osteocalcin
Calcitonin
25-OH Vitamin D
PTH

Cardiac Markers

MYO
cTnl
CK-MB
NT-proBNP
hs-cTnT
BNP
H-FABP*

Endocrine Hormone

LH
FSH
PRL
Testosterone
E2
PRG
DHEA-S
SHBG
17 α -OHP
AMH

Diabetes

Insulin
C-peptide

Growth Hormone

IGF-1
HGH
IGFBP-3

Hepato Fibrosis

PIIINP
Col IV
LN
HA
CG

HBV

HBsAg
Anti-HBs
HBeAg
Anti-HBe
Anti-HBc
HBV PreS1-Ag
Anti-HBc IgM

Rheumatoid Arthritis

Anti-CCP IgG
RF*
RF IgM*
RF IgG*
RF IgA*

Connective Tissue Disease

ANA*
Anti-dsDNA IgG
Anti-SS-B IgG*
Anti-Sm IgG*
Anti-rib P IgG*
Anti-SS-A/Ro IgG*
Anti-Ro-52 IgG*
Anti-Scl-70 IgG*
Anti-Jo-1 IgG*
Anti-nRNP/Sm IgG*

Antiphospholipid syndrome

Anti-β 2-GP1*
Anti-β 2-GP1 IgA*
Anti-β 2-GP1 IgG*
Anti-β 2-GP1 IgM*
Anti-CL*
Anti-CL IgA*
Anti-CL IgG*
Anti-CL IgM*

Autoimmune Liver Disease

Anti-gp210 IgG*
Anti-sp100 IgG*
AMA-M2 IgG*

Autoimmune Diabetes

IAA IgG*
GADA*
IA-2A IgG*
ICA IgG*
ZnT8A IgG*

Inflammation Monitoring

hs-CRP
PCT
IL-6
D-Dimer
Lp-PLA2*

Tuberculosis

TB-IGRA

TORCH

Toxo IgM
CMV IgM
Rubella IgM
HSV-1 IgM
HSV-2 IgM
Toxo IgG
CMV IgG
Rubella IgG
HSV-1 IgG
HSV-2 IgG

Pre-natal

Free β -HCG
P-AFP
β -HCG
Ue3
PAPP-A
Inhibin A*

Hypertension

Aldosterone
ACTH
Cortisol
Renin
Angiotensin II

Allergy

tlgE

Anemia

Folate
Vitamin B12
Ferritin

Thyroid

T3
T4
TSH
FT3
FT4
Anti-TG
Anti-TPO
TG
PTH
rT3
TRAb

*Coming Soon

Autobio Diagnostics Co., Ltd.

No.87, Jingbei Yi Rd, National Eco & Tech Zone, Zheng zhou City, China-450016

Facebook: @autobioGlobal, Youtube Channel: Autobio Diagnostics,

E-mail: info@autobio-diagnostics.com, Whatsapp: +8618595767835, Tel: +86-371-6798-5313

