

SEMI AUTO BIOCHEMISTRY ANALYZER



MODEL: BIOLYSER 4000



MODEL: DI 50

SPECIFICATIONS

Absorbance Range	0.500 – 3.500Abs
Resolution	0.001 Abs (displayed), 0.0001Abs (calculated)
Light Source	Halogen lamp, 12V/20W
Wavelengths	340, 405, 500, 546, 620nm, 3 more ters optional
Wavelength Accuracy	± 1nm
Band Width	< 8nm
Flow Cell	Metal – quartz ow cell
Temperature Control	25°C, 30°C, 37°C ±0.1 and ambient temperature
Carryover	< 1%
Display	LCD display
Interface	RS-232
Output	Builtin Internal and external printer (optional)
Power Supply	AC 100V – 240V ± 10%, 50-60Hz
Net Weight	7Kg
Dimensions L × W × H (mm)	370 x 284 x 318