

ELISA MICROPLATE READER



MODEL: BIOMAX 6500

FEATURES

- Semi automatic
- Easy Windows operation system with touch screen or mouse, large LCD display
- 8-channel optical fiber system enables 5 sec reading for the whole plate
- Bichromatic measurement, calculation modes including ABS, Cut-off, Curve, Linear,
- Log and Exponential regression
- Random positive and negative control setting
- Multi-assay enables up to 12 different assays on one plate
- Large memory, up to 100 test protocols and 10,000 can be stored
- Powerful QC function: Grubs, Westguard Multi-rule, Levey jennings plot (optional)
- Plate shaking
- Auto power-on self-check
- Additional lamp saving
- Multiform result output including patient comprehensive report
- Ability to communicate with PC for data management(optional)
- Easy maintenance