

## ESR Analyzer



**Model: BSR20**

### Specification

- Worker mode : Semi-automatic
- Measuring time : 30 Minutes
- Capacity : 20 pcs
- Maximum measuring speed : 40 pcs/h
- Mechanical/optical precision of detection :  $\pm 0.1\text{mm}$
- Reproducibility of analysis :  $<3\%$  CV
- Redding inaccuracy among each channel :  $<6\%$  CV
- Memory capacity : 264
- Communication interface : USB and optional RS232
- Printer type : Built-in printer
- Automatic temperature conversion :  $18^{\circ}\text{C}$
- Operation part : membrane keyboard operation and barcode scanning port.

